

Date: March 9<sup>th</sup>, 2021

Canadian Heart Rhythm Society  
Device Committee

---

**ADVISORY:**

**Potential premature battery depletion in a subset of  
Biotronik ICD and CRT-D devices (Idova, Iforia, Ilesto / Inventra, Iperia, Itrevia /  
Ilivia, Inlexa, Intica / Ilivia Neo, Intica Neo)**

---

**Nature of the Advisory:**

Premature battery depletion has been seen in a subset of these devices, which have been in distribution since 2013. This results from “lithium plating” which represents a certain mode of lithium deposition on the anodes of the batteries causing a higher rate of battery drain. No serious injury or death have been reported so far.

**Scope of the problem:**

This premature battery depletion could result in loss of high-voltage or pacing therapy. Inability to interrogate or transmit data can also be a manifestation of this problem. The mean interval from ERI to EOS could be reduced to a median of **58 days for high-voltage devices and 6 months for low voltage devices**. The onset for devices to experience this issue is about 2 years with a failure rate estimated at **0.0012%** with a 5-year probability of **0.17%**. The estimated risk for loss of high-voltage therapy is 0.0069% and 0.0015% for loss of pacing capability.

In Canada, there were 24 affected devices on consignment and 388 devices implanted.

### **Response of the CHRS Device Committee:**

- As part of this formal advisory, we recommend that patients be notified about this issue
- All non implanted affected devices should be removed from hospital inventory and returned to Biotronik
- Patients should be followed every 3 months remotely or in person
- No prophylactic replacement is recommended at this stage
- All patients should be enrolled on the Biotronik home monitoring system where possible. This recommendation applies to all patients with CIEDs, but is especially important in this advisory patient group
- In cases where unexpected ERI is reached or inability to interrogate the device or to transmit data is encountered, replacement should be performed:
  - For pacemaker dependant patients or secondary prevention patients, they should be **admitted and replacement performed promptly**
  - In non-pacemaker dependent patients, primary prevention with no previous therapy delivered, replacement should be planned **in ≤ 1 week**.
- BIOTRONIK is working on a software update that reduces the probability of batteries developing this form of lithium plating. Please follow-up with your Biotronik representative for an update on the release date
- The CHRS device committee may update these recommendations should more data become available.

François Philippon, MD, FRCPC, FHRS, FCCS  
Chair, Device Committee

Larry Sterns, MD, FRCPC, FHRS, FCCS  
Past chair, Device Committee

Jason Andrade, MD, FRCPC  
Derek Exner, MD, FRCPC, FHRS  
Clarence Khoo, MD, FRCPC  
Ratika Parkash, MD, FRCPC, FHRS  
Calum Redpath, MBChB, MRCP (UK), PhD  
Raymond Yee, MD, FRCPC, FHRS